# Pocket BVM User Manual

**MANUFACTURED BY:**

Micro BVM™ Systems Ltd.

11 Ruzin St. Jerusalem 93870 Israel

+972-2-6517191

info@microbvm.com

www.microbvm.com

MD

EN

FR

DE

ES

NL

IT

PL

DA

NO

CS

SV

FI

IS

**POCKET BVM USER MANUAL****INTENDED USE****Adult single patient use resuscitator including face mask, reservoir bag and oxygen inlet.**

The **Pocket BVM** Resuscitator is intended for manual pulmonary resuscitation and emergency respiratory support of adult patients.

The Pocket BVM is a single-patient-use device and, as such, no service is required. It should not be re-used on other patients.

**The Pocket BVM is not sterile. Its lifetime is 5 years.**

**The Pocket BVM is to be used by trained personnel only.**

**U.S. Federal law restricts this device to sale by or on the order of a licensed healthcare practitioner.**

**WARNINGS:**

- Oil or grease should not be used in close proximity to oxygen equipment; do not smoke or use open flames when oxygen is in use. Fire may result in both cases.
- Reduced or insufficient airflow may result in brain damage to the patient being ventilated.
- Improper function could lead to inadequate ventilation, resulting in serious injury or death.
- If efficient ventilation cannot be obtained with this device, the operator should immediately switch to mouth-to-mouth and/or mouth-to-nose resuscitation.
- Use in contaminated environments can be hazardous unless entrainment of contaminants is prevented.
- Oxygen flow should not exceed 35 L/min.

**PRECAUTIONS & CAUTIONS:**

- Only persons with adequate training should use the Pocket BVM.
- Users of this product must be familiar with the important information in the directions for use.
- Test each unit prior to use.
- This product is intended for single use only. Do not attempt to clean, decontaminate or sterilize.
- Do not use the resuscitator in toxic or hazardous atmosphere, as silicone rubber is highly permeable to gasses.
- Never store the resuscitator in a squeezed or compressed state, other than folded at the folding zones in the original packaging. Permanent deformation of the bag will occur, which may reduce the ventilation efficiency.
- Avoid contamination with oil and grease. Such products may affect the integrity of the resuscitator materials.

## POCKET BMV USE

1) Unpack and assemble the Pocket BMV from the protective box (Figure 1 below).


Figure 1


### 2) Functional Test for the Pocket BMV, prior to use:

- Check the resuscitator by squeezing a few times to confirm that the Inlet Valve and Patient Valve are both functioning, i.e., that air is entering into the squeeze bag and exiting from the patient valve.
- Check that the bag fills. This is confirmed by visible expansion and relaxation of the bag.
- Using your thumb, block the Patient Port, squeeze the bag and hold for positive pressure, to verify that the system is leak free.
- If there is air leak from the patient valve, tighten the Expiratory connector and try again.
- If any of the components (such as the valves) are not fully functional, do not use this device.

### 3. Ventilate the patient:

- If resuscitating with Oxygen, attach the oxygen supply tube to a regulated oxygen source and set desired flow. Note: O<sub>2</sub> should be administered as soon as possible.
- Use the recommended techniques to position the patient correctly, to open and clear the mouth and airways, and to hold the mask firmly against the face.
- Ensure that the facemask forms an adequate seal with the patient's face. A tight seal between face and mask is essential. Failure to attain a tight seal may result in reduced or no ventilation.
- The Pocket BVM functions must be visible and audible: Listen for expiratory airflow from the valves.
- Ventilate the patient: During insufflation observe the rise of the patient's chest. Release the bag abruptly, listen for the expiratory flow from the patient valve and observe the lowering of the chest. If continued resistance to insufflation is encountered, check for airway obstruction, or correct the backward head tilt.
- If the patient vomits during mask ventilation, immediately clear the patient's airways of vomitus using recommended techniques.
- Clear the Pocket BVM of vomitus, and freely compress the bag a few times before resuming ventilation.
- The correct ventilation frequency may vary. Please follow the current ventilation frequency recommended by regional, national or international guidelines.

### FOLDING AND PACKING FOR STORAGE

**SPECIFICATIONS:**

The body mass range for which the resuscitator is suitable for use	>40 kg (88 lbs.)
Operating environmental limits	-18° C (0° F) to 50° C (122° F) ↓
Storage environmental limits	-40° C (-40° F) to + 60° C (140° F) ↓
Substances, other than air, which can be administered by the resuscitator	Oxygen
Delivered oxygen concentrations under various test conditions	55% (at 2 L/min) to 100% (at 8 L/min)
Gas inlet tube connection	15 mm length x 6 mm O.D.
Stroke-volume range for operator-powered resuscitators	800 ml
Apparatus dead space (backward leakage and forward leakage where appropriate)	<6 ml
Expiratory resistance	3.4 cm H <sub>2</sub> O
Inspiratory resistance	-4.5 cm H <sub>2</sub> O
The end-expiratory pressure generated by the resuscitator in normal use, if greater than 0,2 kPa (2 cmH <sub>2</sub> O)	(<0,2 kPa)
External dimensions of the resuscitator	217 mm (length) x 121 (diameter)
External dimensions of the resuscitator in case	63/73 mm (length) x 73 diameter)
Mass of the resuscitator, face mask and the resuscitator case	Resuscitator 273 g Case. 90 g, Mask 87 g. Total 450 g

**EXPLANATION OF SYMBOLS:**

 Manufacturer	 Year of manufacture of active devices
 European Authorized Representative	 Limit temperature
 Lot number	 Medical device
	 Use by

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